


MicroHub® is almost here!

The modern Lab Intelligence Platform

Finally manage the lab, operate your routine, and empower data from **one single software**.

AI powered **validation**, customizable **dashboards, statistics**, and an **integration-first** design are just some of the features you will love.



Last phases before launch

MicroHub® is almost ready, get in contact in advance!

MicroHub® is about to be released on the market. **Forward thinkers** know that this is the time to get in touch and start to **recover informations** for future implementations. Browse this brochure or even better, get in touch!

Request a demo!

No brochure will tell you everything you need to know: **get in touch with us**, we love to talk to people!

Write to:

sales.wasp@copangroup.com

Or visit our website...

www.copangroup.com →

...and click

Contact us



MicroHub®

Briefly said, MicroHub® is a gamechanger

There is no plain definition to what MicroHub® is: it connects all the lab automations and softwares **like a middleware**, looks like an L.I.S. as it works as an interface with patient, sample, test and workup data, it aggregates and analyzes them like a **management platform**, and finally it streamlines the validation process like... no other, actually.

Please welcome the real **Microbiology Hub**.

Why is it unique?

Connection without compromises

MicroHub® is not the usual "We connect our own automation" software. It is designed to jump **right in the middle** of your microbiology lab automation environment, bringing every available data together in real-time right into your hands. **One single interface** to work with L.I.S. as well as all platforms, workloads, and results, manage the validation phase in a flash.



Leverage MicroHub® to...




Advance your report phase to the next level with serial-, mass-, and auto-validation!



Collect and arrange all the data the way you really need.



Monitor all lab automation status and workload in one single view.



Analyze all the lab activities with statistics and aggregated data.

No time to read?

Jump right to the second half of this brochure to peek into screenshots and feature highlights.



What will MicroHub® do for you?


Connect and collect to a common base

Main priority: retrieve all the data

Right now in your lab, patients, samples, tests, and automation are generating a cloud of data that awaits to express its potential. MicroHub® integrates with **L.I.S** and all the other **digitally connected automated systems**, and can gather information ranging from patient anagraphics, to test results, plate and tubes id, automation processing status, and obviously intermediate outcomes. Basically anything that gets recorded.




A real-time Data Stack



Keystones



*Multichannel L.I.S. HL7
Connection*



*Different inputs brought
to a common base*



*Modulate your UI over the info
you want to get together*

Make them accessible

Your own priority at first-sight

Imagine every information to be connected in a huge shapeless cloud, or less dramatically, in an impossibly complex spreadsheet. The more data you bring in, the more the spreadsheet grows. How does MicroHub® enable its readability?

Flexible User Interface, widgets and dashboards, object filtering.

Just focus on your interpretation key, and get your data arranged by **your point of view**.

Patient First?

Sort all samples and tests for each patient's status, historical clinical data, or outcome. Perfect for **cross-talks** between departments when **facing specific clinical cases**.

Sample First?


Breakdown any sample downstream items (plates, ID targets, secondary tubes, enrichment tubes, and so on) and all the tests performed with it, their status and their results.

Management First?

Epidemiology, workload, quality, TAT. Perfect for **reporting and analyzing** overall lab performances and planning for future **improvements**.

Quality First?

Explore every step of your lab workflows, check for **inconsistencies, errors** and address **black spots**, for QA-proof lab escalation.




Do-it-yourself customization!

You read it half a dozen times:

MicroHub® is a software designed to address the needs of **many different professionals**. That's why a **self-customization** is a feature that cannot be missed here.

- Create your own widgets for easy access to your activity lists and share it with your colleagues
- Organize your dashboard to highlight priorities or generic topics
- Manage different editor rights
- Append notes so everyone knows exactly the purpose of each widget

All of this without depending on your provider service!



We love algorithms

A new way of thinking **Technical** and **Clinical Validation**

Like our top-notch softwares PhenoMATRIX® and Radian® Expert System, MicroHub® leverages Artificial Intelligence principles to fight **against repetitive and error-prone tasks**. Build your rule set, enjoy MicroHub® bringing up all the tests that match those rules, and then... choose the way you prefer to proceed.



Serial Validate

In a single page all the info you need to take action: one click to report directly to L.I.S. and pass to the next page.




Mass-validate

Just one click... and every test that complies with your validation rules is sent to L.I.S.



Auto-validate

Forget about it. Everything that fits your rules gets validated right away.



L.I.S.

We love algorithms even more

A new way of thinking **Workup Waterfall**

On the same principle, the use of expert rules can be applied not just for validation! Thanks to the wide range of data retrieved, MicroHub® can also automatize **chunks of the lab workflow**, especially the established ones that imply multiple lesser steps.

A simple example for a complex concept (tongue twister included)

Positive Blood Culture to Gram Slide

Ouch! We have a positive blood culture! A pretty established practice is to immediately go for a morphological and staining examination for the patient positive blood cultures with a slide. Why then a lab tech should check for the positivity, get the tube, scan the label, submit an order (maybe print some document) and finally perform it... instead of just doing it?

With a simple rule set, MicroHub® will check the positivity, and append your test of choice to the sample! Total time required? Few seconds in the worst case.



Getting the maximum from WASPLab®

WASPLab® and MicroHub® speak WISELY

Up to now we said that MicroHub® is able to get the best from all the lab automations. It's obvious thought that it is a Copan Software designed to **work in perfect synergy** with Copan Automation, as a perfect expression of our WISE approach.

Expanded control over WASPLab® orders

Enhanced WASPLab® workup Automatization

Optimal Radian® Expert System integration

Request a demo!

No brochure will tell you everything you need to know: **get in touch with us**, we love to talk to people!

Write to:

sales.wasp@copangroup.com

Or visit our website...

www.copangroup.com →

...and click

Contact us

MicroHub[®] hands-on

Let's peek into some screenshots

The core of patient-centric features

This is the one page you'll use to investigate all the sample/test/workup/result layers **connected to a patient**. On the left you will find patient ID and demographics, you can choose the samples and related tests for each.

In the middle all the core information and on the right all the horizontal ones, that you want to have right at your fingertips to take the right decision.

Ready to proceed? Validate the test or just add comments and send the result to L.I.S.

Unique page for all the patient information

The screenshot displays the Microbiology Validation interface for patient **MARSHALL, Bill** (ID: IDP0071). The interface is divided into several sections:

- Left Sidebar:** Contains navigation icons for ACTIVITIES, VALIDATION, MONITOR, STATISTICS, RELEASE, and IDENTIFICATIONS.
- Top Header:** Shows the patient's name, ID, and a list of activities (MICROBIOLOGY VALIDATION).
- Central Panel:** Displays test results for **MICETI - YEASTS CULTURE**. It includes a table of isolates with columns for CODE, TI, CREATION DA., TI, MICROORGANISM, BACTERIAL, WORKUPS, ACTIVITY ST., and ACTIONS. Two isolates are listed: *escocci - Escherichia coli* and *bianche - Bianche*.
- Right Panel:** Contains a list of other information, including ALARMS, CLINICAL INFO, SAMPLE INFO, ORDER DETAIL, BOTTLES, MEDIA, and OTHER LIS RESULTS.

Two callouts highlight specific features:

- Attribution of isolate type: pathogen/contaminant** (pointing to the isolate table).
- Parallel comparison between tests from different samples** (pointing to the test results section).

Clinical evaluation: full sample, individual test, validation

The screenshot displays the Microbiology Validation interface for patient **MARSHALL, Bill** (ID: IDP0071). The interface is divided into several sections:

- Left Sidebar:** Contains navigation icons for ACTIVITIES, VALIDATION, MONITOR, STATISTICS, RELEASE, and IDENTIFICATIONS.
- Top Header:** Shows the patient's name, ID, and a list of activities (MICROBIOLOGY VALIDATION).
- Central Panel:** Displays test results for **URC - URINE CULTURE**. It includes a table of isolates with columns for CODE, TI, CREATION DA., TI, MICROORGANISM, BACTERIAL, WORKUPS, ACTIVITY ST., and ACTIONS. One isolate is listed: *escocci - Escherichia coli*.
- Right Panel:** Contains a list of other information, including ALARMS, CLINICAL INFO, SAMPLE INFO, ORDER DETAIL, BOTTLES, MEDIA, and OTHER LIS RESULTS.

Two callouts highlight specific features:

- Configurable alert messages** (pointing to the ALARMS section).
- Labtech operations: test management, workups, comments** (pointing to the isolate table).

Technical and Clinical Worklists, all configurable!

To understand real the power of this tool you should not just think about the sample type though: you can build worklist based on **test type, workup type, workup status, sample status, and so on**.


Don't limit your requirements on technical informations, go also for the clinical ones, and don't forget that you can customize your whole dashboard, not just single widgets!



Lab Monitor Dashboard

Check the real-time monitor and manage!

Recall in a single dashboard the workload of blood cultures, streaking and incubation, IDs and ASTs workups, and then dig deeper to see graphics and tables to get a new perspective on your lab daily routine.



Get the best from ID integration

Identification earns a dedicated page. From here you will have the full **breakdown of isolate identification data**, ranging from material isolate number, presumptive result, and the quality rate of the identification.

On the right are the actions you may want to perform: from your own workstation you can confirm, forward to LIS, comment... or even ask for a repeat or new workup.

Single or double spot management with extended information

Selective results sending to LIS (with preselected results)

Management of raw Maldi Results (Bruker)

Manage re-runs and confirmed results

PATIENT	SAMPLE	MATERIALE	ISOLATE	PRESUMPT	SPOT	MICROORGANISM	QUALITY	CONFIRM	SEND TO LIS	REPEAT	IGNORE	COMMENTS
BLESCE, ANTONIO	B2022092701	UR - Urine	1		A2	klepne - Klebsiella pneumoniae	●	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
NUZZO, BIAGIO	B2022092702	UR - Urine	1		A3	klepne - Klebsiella pneumoniae	●	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
SCHIRRU, DANIELE	B2022092703	UR - Urine	1		B1	staeppi - Staphylococcus epidermidis	●	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
GRISOLIA, FILOMENA	B2022092704	UR - Urine	1		B2	klepne - Klebsiella pneumoniae	●	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
RICCA, ANGELO	B2022092705	UR - Urine	1		B3	esccol - Escherichia coli	●	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
RATTACASO, ANGELA	B2022092706	UR - Urine	1		B4	esccol - Escherichia coli	●	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
BOTTA, GIORGIA	B2022092708	UR - Urine	1		C3	esccol - Escherichia coli	●	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
MOMBELLI, VIRGINIA	B2022092709	UR - Urine	1		C4	entclo - Enterobacter cloacae	●	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Consult both evaluated or to-be-evaluated lists

PROBABILISTIC RESULTS

CONFIRM SELECTED

PROGRAMMING DATE/TIME:

RESULTS (8)

PATIENT	SAMPLE	MATERIALE	ISOLATE	PRESUMPT	SPOT	MICROORGANISM	QUALITY
BLESCE, ANTONIO	B2022092701	UR - Urine	1		A2	klepne - Klebsiella pneumoniae	MEDIUM
NUZZO, BIAGIO	B2022092702	UR - Urine	1		A3	klepne - Klebsiella pneumoniae	HIGH
SCHIRRU, DANIELE	B2022092703	UR - Urine	1		B1	staeppi - Staphylococcus epidermidis	
GRISOLIA, FILOMENA	B2022092704	UR - Urine	1		B2	klepne - Klebsiella pneumoniae	
RICCA, ANGELO	B2022092705	UR - Urine	1		B3	esccol - Escherichia coli	
RATTACASO, ANGELA	B2022092706	UR - Urine	1		B4	esccol - Escherichia coli	
BOTTA, GIORGIA	B2022092708	UR - Urine	1		C3	esccol - Escherichia coli	
MOMBELLI, VIRGINIA	B2022092709	UR - Urine	1		C4	entclo - Enterobacter cloacae	

RESULTS (10)

MICROORGANISM	SCORE	QUALITATIVE PROBABILITY
klepne - Klebsiella pneumoniae	1.32	+
klepne - Klebsiella pneumoniae	1.31	+
klepne - Klebsiella pneumoniae	1.29	+
klepne - Klebsiella pneumoniae	1.27	+
klepne - Klebsiella pneumoniae	1.21	+
klepne - Klebsiella pneumoniae	1.13	+
kleoxy - Klebsiella oxytoca	1.07	-
kleoxy - Klebsiella oxytoca	1.03	-
kleoxy - Klebsiella oxytoca	1.02	-
kleoxy - Klebsiella oxytoca	1.02	-

RESULTS (10)

MICROORGANISM	SCORE	QUALITATIVE PROBABILITY
klepne - Klebsiella pneumoniae	2.67	+++
klepne - Klebsiella pneumoniae	1.65	++
klepne - Klebsiella pneumoniae	1.47	++
klepne - Klebsiella pneumoniae	1.32	+
klepne - Klebsiella pneumoniae	1.31	+
klepne - Klebsiella pneumoniae	1.23	+
kleoxy - Klebsiella oxytoca	1.09	-
kleoxy - Klebsiella oxytoca	1.07	-
kleoxy - Klebsiella oxytoca	1.02	-
kleoxy - Klebsiella oxytoca	1.02	-




Copan

Not just here-and-now

Often underestimated, the **retrospective analysis of data** brings a huge added value and a solid base to build up your lab proficiency and impact.

If you are managing the lab at an higher level, peek into the **Statistics panel**.

Explore epidemiology data, monthly based workloads, turnaround times, overview process statistics, or even breakdown your Quality KPI, like contamination rates or microbiology tests overconsumption.



Workflow Integrated System Environment

We fell in love with Microbiology a long time ago.

Since then we worked tirelessly with one goal in mind: designing preanalytics as a whole.

Designing products that empower each other to address the most critical phase^[1] of the microbiological analysis. Every product in the Copan range, finds its reason-to-be in this wider wiser vision.



[1] Lippi Giuseppe, Chance Jeffrey J., Church Stephen, et al. "Preanalytical quality improvement: from dream to reality" Clinical Chemistry and Laboratory Medicine, vol. 49, no. 7, 2011, pp. 1113-1126.



This document may contain product information otherwise not accessible or valid in your country. Please be aware that Copan Wasp s.r.l. does take any responsibility for accessing such information which may not comply with any valid legal process, regulation, registration or usage in the country of your origin. Product clearance and availability restrictions may apply in some Countries. Please refer to Copan website (www.copangroup.com) to view and/or download the most recent version of the brochure. This document is mainly intended for marketing purposes, always consult product insert for complete information. The use of this product in association with diagnostic kits or instrumentation should be internally validated by the user. ©2022 Copan Italia. All rights reserved. All product names, trademarks and registered trademarks mentioned herein are property of their respective owners.
Code: MMWB-MHBooRoo 2022-12



@copangroup

Copan Wasp s.r.l.

Via A. Grandi 32,
25125 Brescia, Italy

t | f +030 2687218

@ | wasp@copangroup.com

www.copangroup.com