

# PASTEUR PIPETS

SAFE TRANSFERRING AND DISPENSING OF DIFFERENT KIND OF LIQUIDS



## SINGLE-USE LABORATORY PRODUCTS

Disposable plastic Pasteur pipets, commonly referred to as Transfer Pipets, are essential products widely used by hospital, research and industrial laboratories for transferring liquids and clinical samples during diagnostic testing or analytical research. COPAN Disposable Transfer Pipets are ideal for easily transferring and dispensing liquids\* in all types of laboratory.

- Ideal for **transferring** and **dispensing liquids** in all types of laboratories:
  - blood culture
  - urine
  - bacteriology culture
  - dairy products
  - different types of viscose liquids
- **Sterile** pipets available in either single or multiple peel-pouch formats
- Available in a wide range of standard and designed types

## COPAN - PASTEUR PIPETS



Total Length	Tip Diameter	Drops/ml	Total Capacity	Graduated up to	Grad. Intvl
150mm	3.3mm	21	7ml	3ml	0.5ml

Cat.No	Description	Pack Size
200C	non sterile graduated pipet	6x500pcs per case
200CS01	sterile graduated pipet, individually wrapped	6x500pcs per case
200CS05	sterile graduated pipet, 5 pieces per pack	5x100pcs per case
200CS10	sterile graduated pipet, 10 pieces per pack	10x50pcs per case



Total Length	Tip Diameter	Drops/ml	Total Capacity	Graduated up to	Grad. Intvl
151mm	3.0mm	24	5ml	1ml	0.5ml

Cat.No	Description	Pack Size
201C	non sterile graduated pipet	6x500pcs per case
201CS01	sterile graduated pipet, individually wrapped	6x500pcs per case
201CS05	sterile graduated pipet, 5 pieces per pack	5x100pcs per case
201CS10	sterile graduated pipet, 10 pieces per pack	10x50pcs per case



Total Length	Tip Diameter	Drops/ml	Total Capacity	Graduated up to	Grad. Intvl
227mm	3.0mm	25	6ml	1.5ml	0.25ml

Cat.No	Description	Pack Size
204C	non sterile graduated pipet	6x500pcs per case
204CS01	sterile graduated pipet, individually wrapped	6x500pcs per case
204CS05	sterile graduated pipet, 5 pieces per pack	5x100pcs per case
204CS10	sterile graduated pipet, 10 pieces per pack	10x50pcs per case



Total Length	Tip Diameter	Drops/ml	Total Capacity	Graduated up to	Grad. Intvl
136mm	2.5mm	25	3ml	1ml	0.25ml

Cat.No	Description	Pack Size
218C	non sterile graduated pipet	6x500pcs per case

## COPAN - PASTEUR PIPETS



Total Length	Tip Diameter	Drops/ml	Total Capacity
82mm	3.0mm	24	4ml

Cat.No	Description	Pack Size
203C	non sterile narrow stem pipet	1000pcs per box



Total Length	Tip Diameter	Drops/ml	Total Capacity
147mm	2.3mm	28	4.4ml

Cat.No	Description	Pack Size
202C	non sterile narrow stem pipet	6x500pcs per case
202CS01	sterile narrow stem pipet, individually wrapped	6x500pcs per case
202CS05	sterile narrow stem pipet, 5 pieces per pack	5x100pcs per case
202CS10	sterile narrow stem pipet, 10 pieces per pack	10x50pcs per case



Total Length	Tip Diameter	Drops/ml	Total Capacity
61mm	2.5mm	27	1.4ml

Cat.No	Description	Pack Size
209C	non sterile narrow stem pipet	6x500pcs per case



Total Length	Tip Diameter	Drops/ml	Total Capacity	Graduated up to	Grad. Intvl
153mm	2.7mm	24	6ml	2ml	0.5ml

Cat.No	Description	Pack Size
212C	non sterile blood bank pipet	10x500pcs per case



Total Length	Tip Diameter	Drops/ml	Total Capacity
103mm	1.0mm	54	1.6ml

Cat.No	Description	Pack Size
222C	non sterile fine tip pipet	6x400pcs per case
222CS01	sterile fine tip pipet, individually wrapped	6x400pcs per case
222CS10	sterile fine tip pipet, 10 pieces per pack	10x50pcs per case



Total Length	Tip Diameter	Drops/ml	Total Capacity
147mm	1.5mm	48	9ml

Cat.No	Description	Pack Size
223C	non sterile fine tip pipet	6x400pcs per case
223CS01	sterile fine tip pipet, individually wrapped	6x500pcs per case
223CS05	sterile fine tip pipet, 5 pieces per pack	5x100pcs per case
223CS10	sterile fine tip pipet, 10 pieces per pack	10x50pcs per case



Total Length	Tip Diameter	Drops/ml	Total Capacity
147mm	1.0mm	58	5.8ml

Cat.No	Description	Pack Size
224C	non sterile fine tip pipet	6x500pcs per case
224CS01	sterile fine tip pipet, individually wrapped	6x500pcs per case
224CS10	sterile fine tip pipet, 10 pieces per pack	10x50pcs per case



Total Length	Tip Diameter	Total Drops	Capacity	Bulb Draw
131mm	1.5mm	10	0.7ml	0.2ml

Cat.No	Description	Pack Size
220C	non sterile paddle pipet	6x500pcs per case

## COPAN - PASTEUR PIPETS



Total Length	Tip Diameter	Drops/ml	Total Capacity
86mm	5.0mm	25	1.5ml
Cat.No	Description	Pack Size	
205C	non sterile general purpose pipet	6x1000pcs per case	



Total Length	Tip Diameter	Drops/ml	Total Capacity
123mm	2.6mm	24	2.3ml
Cat.No	Description	Pack Size	
211C	non sterile general purpose pipet	10x600pcs per case	



Total Length	Tip Diameter	Drops/ml	Total Capacity
142mm	2.8mm	26	4.5ml
Cat.No	Description	Pack Size	
219C	non sterile general purpose pipet	6x500pcs per case	




Total Length	Tip Diameter	Drops/ml	Total Capacity
167mm	4.0mm	17	15ml
Cat.No	Description	Pack Size	
207C	non sterile extra large pipet	1200 pieces	



Total Length	Tip Diameter	Drops/ml	Total Capacity
290mm	3.0mm	21	24ml
Cat.No	Description	Pack Size	
227C	non sterile extra large pipet	6x100pcs per case	

# DROPLESS PIPET, COUNTING DROP NO MORE!




> Simplified rapid membrane test protocol

> Accurate one-step volume delivery

## SIMPLY FILL THE PIPET AND RELEASE!

**FILL:** Squeeze and release bulb to draw up liquid. Excessed liquid is siphoned automatically in two overflow tanks.

**RELEASE:** Squeeze bulb another time to release precise volume just off the stem of the pipet. Proceed until the stem is fully empty.



Cat.No	Description	Pack Size
214C	100ul pipet	2x500pcs in 6 boxes

\*The use with solvents (polar, non-polar, ionic) or corrosive/water-reactive liquids shall be validated by the user.

**Copan**  
innovating together

This document may contain product information otherwise not accessible or valid in your country. Please be aware that Copan Italia S.p.A. does take any responsibility for accessing such information which may not comply with any valid legal process, regulation, registration or usage in the country of your origin. Product clearance and availability restrictions may apply in some Countries. Please refer to Copan website ([www.copangroup.com](http://www.copangroup.com)) to view and/or download the most recent version of the brochure. This document is mainly intended for marketing purposes, always consult product insert for complete information. The use of these products in association with diagnostic kits or instrumentation should be previously validated by the user. ©2021 Copan Italia. All rights reserved. The trademarks mentioned herein are property of Copan Italia S.p.A.  
Code: JMKCo25R01



@copangroup

**Copan Italia s.p.a.**  
Via Francesco Perotti 10,  
25125 Brescia, Italy

t | f +030 2687211  
@ | [info@copangroup.com](mailto:info@copangroup.com)  
[www.copangroup.com](http://www.copangroup.com)