

FLOQSwabs®

IN PEELPOUCH



Code	Description	Break Point ¹	Specimen ²
519CS01	REGULAR 	100 mm	Nasal, throat, vaginal, groin, armpit, rectal, wound, buccal, faeces
502CS01	REGULAR 	80 mm	Nasal, throat, vaginal, groin, armpit, rectal, wound, buccal, faeces
520CS01	REGULAR 	30 mm	Nasal, throat, groin, armpit, wound, buccal, faeces
5U060S01	REGULAR 	no breaking point	Nasal, throat, vaginal, groin, armpit, rectal, wound, buccal, faeces
537CS01	REGULAR, CONIC TIP 	no breaking point	Nasal, throat, vaginal, groin, armpit, rectal, wound, buccal, faeces
518CS01	MINITIP 	100 mm	Urethral, nasal, pediatric
501CS01	MINITIP 	80 mm	Urethral, nasal, pediatric

FLOQSwabs® - STERILE | Packaging: 1000 pcs (10 boxes of 100 pcs in a case)



Code	Description	Break Point ¹	Specimen ²
525CS01	ULTRA-THIN MINITIP 	80 mm	Urethral
503CS01	FLEXIBLE 	100 mm	Naso-pharyngeal
5U074S01****	FLEXIBLE, IN BULK 	100 mm	Naso-pharyngeal
534CS01	FLEXIBLE 	no breaking point	Naso-pharyngeal
5U073S01****	FLEXIBLE, IN BULK 	no breaking point	Naso-pharyngeal
516CS01	MINI FLEXIBLE 	100 mm	Neonatal
535CS01	MINI FLEXIBLE 	no breaking point	Neonatal

FLOQSwabs® - STERILE | Packaging: 1000 pcs (10 boxes of 100 pcs in a case) | **** Packaging: 1000 pcs in a case

FLOQSwabs®

IN PEELPOUCH




Code	Description	Break Point ¹	Specimen ²
511CS01*	CONE-SHAPED 	80 mm	Endo-esocervical
570CS01*	L-SHAPED 	80 mm	Endo-esocervical
5U053S01*	FINE TIP 	80 mm	Endocervical, faeces
56380CS01	LONG CONTOURED WITH STOPPER 	80 mm	Mid Turbinate (adult)
56780CS01	CONTOURED WITH LOW STOPPER 	80 mm	Mid Turbinate (pediatric)
56750CS01	SHORT CONTOURED WITH STOPPER 	50 mm	Mid Turbinate (pediatric)
5U069S01.20**	FLEXIBLE, BUNDLES OF 20 	no breaking point	Naso-pharyngeal

FLOQSwabs® - STERILE | Packaging: 1000 pcs (10 boxes of 100 pcs in a case) | *Packaging: 500 pcs (10 boxes of 50 pcs in a case) | **Packaging: 1500 pcs in a case

FLOQSwabs®
IN PEELPOUCH










Code	Description	Break Point ¹	Specimen ²
5U068S01.25**	FLEXIBLE, BUNDLES OF 25 	no breaking point	Naso-pharyngeal

FLOQSwabs® - STERILE | **Packaging: 1500 pcs in a case

FLOQSwabs®

IN TUBE










Code	Description	Break Point ¹	Specimen ²
519C	REGULAR 	100 mm	Nasal, throat, vaginal, groin, armpit, rectal, wound, buccal, faeces
552C	REGULAR 	80 mm	Nasal, throat, vaginal, groin, armpit, rectal, wound, buccal, faeces
552C.80*	REGULAR 	80 mm	Self-Vaginal
552C.80PB*	REGULAR, PEELABLE BARCODE 	80 mm	Self-Vaginal
520C	REGULAR 	30 mm	Nasal, throat, groin, armpit, wound, buccal, faeces
5U032S	REGULAR 	no breaking point	Nasal, throat, vaginal, groin, armpit, rectal, wound, buccal, faeces
5U085S	REGULAR 	no breaking point	Nasal, throat, vaginal, groin, armpit, rectal, wound, buccal, faeces

FLOQSwabs® - STERILE | Packaging: 1000 pcs (10 boxes of 100 pcs in a case) | *Packaging: 600 pcs (6 boxes of 100 pcs in a case)

FLOQSwabs®








IN TUBE



Code	Description	Break Point ¹	Specimen ²
5U041S	REGULAR, CONIC TIP 	no breaking point	Nasal, throat, vaginal, groin, armpit, rectal, wound, buccal, faeces
5U051S	REGULAR PLUS, ROUNDED TIP 	no breaking point	Vaginal
5E046S*	REGULAR PLUS, ROUNDED TIP, PEELABLE BARCODE 	no breaking point	Self-Vaginal
5U005S***	REGULAR, DUAL 	no breaking point	Nasal, throat, vaginal, groin, armpit, rectal, wound, buccal, faeces
5U002S	REGULAR WITH LOW STOPPER 	80 mm	Rectal, entero-pathogens
528C	LARGE TIP WITH STOPPER 	no breaking point	Buccal
52980C	CONE-SHAPED TIP 	80 mm	Endo-esocervical, vaginal, others.

FLOQSwabs® - STERILE | Packaging: 1000 pcs (10 boxes of 100 pcs in a case) | *Packaging: 500 pcs (5 boxes of 100 pcs in a case) | ***Packaging: 600 pcs (6 boxes of 100 pcs in a case)



Code	Description	Break Point ¹	Specimen ²
518C	MINITIP 	100 mm	Urethral, nasal, pediatric
551C	MINITIP 	80 mm	Urethral, nasal, pediatric
553C	FLEXIBLE MINITIP 	100 mm	Naso-pharyngeal
534C	FLEXIBLE MINITIP 	no breaking point	Naso-pharyngeal
516C	MINI FLEXIBLE 	100 mm	Neonatal
535C	MINI FLEXIBLE 	no breaking point	Neonatal
560C	ULTRA-THIN MINITIP 	80 mm	Urethral

FLOQSwabs® - STERILE | Packaging: 1000 pcs (10 boxes of 100 pcs in a case)



This document may contain product information otherwise not accessible or valid in your country. Please be aware that Copan Italia S.p.A. does take any responsibility for accessing such information which may not comply with any valid legal process, regulation, registration or usage in the country of your origin. Product clearance and availability restrictions may apply in some Countries. Please refer to Copan website (www.copangroup.com) to view and/or download the most recent version of the brochure. This document is mainly intended for marketing purposes, always consult product insert for complete information. The use of these products in association with diagnostic kits or instrumentation should be previously validated by the user. ©2021 Copan Italia. All rights reserved. The trademarks mentioned herein are property of Copan Italia S.p.A.
Code: JMK003R07