

# Universal Transport Medium (UTM-RT®)

*With new tube format & formal molecular claim*

Copan **Universal Transport Medium (UTM-RT®)** System is intended for the collection and transport of clinical specimen containing viruses, chlamydiae, mycoplasma or ureaplasma from the collection site to the testing laboratory. UTM-RT® can be processed using standard clinical laboratory operating procedures for viral, chlamydial, mycoplasma and ureaplasma culture. UTM-RT® is intended for the stabilization and transportation of an unprocessed upper respiratory clinical specimen suspected of containing respiratory viruses' nucleic acids. **UTM-RT® is intended for use with compatible molecular assays.**

A new format of the well-known UTM-RT® is available. This new product configuration allowed the **direct loading** on Roche cobas® 5800/6800/8800 systems, **supporting the fully validated workflow**. In addition, this UTM-RT® include a formal claim for its use for the molecular detection of upper respiratory pathogens.

## Features of the new format:




- Same UTM-RT® medium of current formats from Copan - no need to medium change management from actual Copan UTM-RT® formats
- 3 ml volume / no glass beads
- Formal molecular claim: UTM-RT® is intended for use with compatible molecular assays.
- CE-IVDR & FDA Cleared
- Suitable for direct loading on the Roche cobas® 5800/6800/8800 platforms
- Designed for Roche cobas 5800/6800/8800, but even compatible and validated with all molecular platforms






## Need more info?

Visit our web-site: <https://www.copangroup.com/>  
or contact us at: [info@copangroup.com](mailto:info@copangroup.com).

**Copan**  
innovating together

Cat N.	Description	Pack size	Sample*
<b>3U090N</b>	15.5x95 mm round bottom tube filled with 3ml UTM® medium - no beads - shelf life: 18 months 	300 pieces (6 boxes of 50 pieces)	N/A
<b>3U091N01</b>	15.5x95 mm round bottom tube filled with 3ml UTM® medium - no beads + Regular FLOQSwabs® - shelf life: 18 months 	300 pieces (6 boxes of 50 pieces)	Regular shape. Possible selection of collection sites: Nose, oropharynx, skin, vagina, rectum
<b>3U092N01</b>	15.5x95 mm round bottom tube filled with 3ml UTM® medium - no beads + Progressive minitip FLOQSwabs® - shelf life: 18 months 	300 pieces (6 boxes of 50 pieces)	Progressive minitip shape. Possible selection of collection sites: Nasopharynx, eye, hear

**Product code available in US:**

Cat N.	Description	Pack size	Sample*
<b>3C098N</b>	15.5x95 mm round bottom tube filled with 3ml UTM® medium - no beads - shelf life: 18 months 	300 pieces (6 boxes of 50 pieces)	N/A
<b>3C099N01</b>	15.5x95 mm round bottom tube filled with 3ml UTM® medium - no beads + Regular FLOQSwabs® - shelf life: 18 months 	300 pieces (6 boxes of 50 pieces)	Regular shape. Possible selection of collection sites: Nose, oropharynx, skin, vagina, rectum
<b>3C100N01</b>	15.5x95 mm round bottom tube filled with 3ml UTM® medium - no beads + Progressive minitip FLOQSwabs® - shelf life: 18 months 	300 pieces (6 boxes of 50 pieces)	Progressive minitip shape. Possible selection of collection sites: Nasopharynx, eye, hear

\*Suggested table. Please refer to your GLP procedures to choose the most appropriate device for the specific sampling site



@copangroup

**Copan Italia s.p.a.**

Via Francesco Perotti 10,  
25125 Brescia, Italy

t | f +030 2687211

@ | info@copangroup.com

www.copangroup.com